



SMARTSENSOR
AR300+
-32 to 400C



HF120 Human Body Infrared Thermometer

- **Product No.:** HF120
-

- SPECIFICATIONS

Measuring range: -32 to 400°C D:S=12:1 Emissivity: 0.1 to 1.00 adjustable Accuracy: $\pm 1.5\%$ or $\pm 1.5^\circ\text{C}$ °C/F Selection: Yes Response time: 500ms Resolution: 0.1°C Data hold function: Yes	
--	--

Laser Target pointer function: Yes Auto Power off: Yes Display: Green Backlight	