

AM5218MZT



The Dino-Lite Edge DVI series AM5218MZT offers stunning image quality with minimal display lag within a wide range of magnification. Connecting directly to a LCD monitor with 720p, the AM5218MZT can be well-suited to a host of inspection or assembly applications where space is a luxury or security is a concern.

## Overview



DVI interface

The DVI interface allows Dino-Lite to deliver real-time images to a LCD monitor without the presence of a computer.



High optical resolution

The superior optics adopted in the Edge series reveals the finest details, answering the needs of the most demanding microscopy applications.



Microtouch button

This touch-sensitive button allows to freeze/unfreeze image or to control the LED on/off with a long-tap.



Scroll Lock

The scroll lock ensures the focus knob staying at the desired focus or magnification position without worry of unintentional movement.



nterchangeable caps

The interchangeable caps provide adaptability to numerous applications with alternative lighting or object interface, such as but not limited to diffused-light, ring-light, coaxial-light, etc.

## Interchangeable front caps



N3C-C / Close Cap
This cap protects the lens and LED lights from contamination of dust, debris,



N3C-D / Diffuser Cap This cap diffuses the LED light.



N3C-D2 / Opal Diffuser Cap This cap diffuses the LED light.



N3C-E / Extended Open Cap Dino-Lite Edge (stand type) will focus at approximately 200x when the cap touches surface.



magnification.



N3C-O / Open Cap
This is the standard cap for normal usage.



N3C-S / Sidelight Cap
This cap creates images with more depth and texture.

Informa	ation about	working dist	ance and f	ield of view
M	WD	FOV (x)	FOV (y)	DOF
20	60.2	19.5	11.7	2.5
30	33.5	13	7.8	1.8
40	20.9	9.8	5.9	1.5
50	13.9	7.8	4.7	
60	9.7	6.5	3.9	
70	7.1	5.6	3.4	1.0
80	5.5	4.9	3	
90	4.5	4.3	2.6	
100	4.1	3.9	2.3	
110	4	3.6	2.1	
120	4.1	3.3	2	
130	4.5	3	1.8	
140	5	2.8	1.7	
150	5.6	2.6	1.6	
160	6.3	2.4	1.5	
170	7.1	2.3	1.4	
180	8	2.2	1.3	
190	8.9	2.1	1.2	
200	9.9	2	1.2	
210	10.9	1.9	1.1	
220	11.9	1.8	1.1	0.1

M = magnification rate WD = working distance (without front cap) FOV = field of view DOF= depth of field Unit = mm

## Specification

Model: AM5218MZT Dino-Lite Edge

Interface : DVI

Product Resolution : HD 720p Magnification : 20x~220x

Frame Rate: 60fps Sensor: Color CMOS Microtouch: Yes

Unit Dimension: 10.5cm (H) x 3.2cm (D)

Unit Weight: 125g