

KERN 326-13

Single weight, ECO-Shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	80 x 80 x 144 mm
Gross weight:	5 kg
Net weight:	5 kg
Shipping method:	Parcel service

Construction

Dimensions:	80 144
Material:	polished stainless steel

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	5 kg
OIML class:	F1
Tolerance:	25 mg

Accessories

Aluminium box:	317-130-600
Plastic box:	347-130-400
Wooden box:	317-130-100

Service

DAkkS Certificate:	962-443
Positive tolerance:	921-443

Category

Brand:	KERN
--------	------