

PCE Americas Inc. 711 Commerce Way Suite 8 Jupiter FL-33458 USA From outside US: +1 Tel: (561) 320-9162 Fax: (561) 320-9176 info@pce-americas.com PCE Instruments UK Ltd. Units 12/13 Southpoint Business Park Ensign way Hampshire / Southampton United Kingdom, SO31 4RF From outside UK: +44 Tel: (0) 2380 98703 0 Fax: (0) 2380 98703 9 info@pce-instruments.com

www.pce-instruments.com/english www.pce-instruments.com

## School Environmental Meter PCE-PH 15 Simple pH-Meter / Resolution 0.01 pH / Removable electrode / Easy operation / Small and handy meter / Onepoint adjustment / For simple measuring tasks

The school environmental meter PCE-PH 15 is a measuring instrument for the determination of the pH value in liquids. Thanks to its ease of use, the PCE-PH 15 school environmental meter is ideal for use in, for example, schools or educational institutions. In addition to the normal measuring function, the school environmental meter has a "HOLD" function. Here it is possible to freeze the currently displayed measured value on the school environmental meter.

With an accuracy of  $\pm$  0.1 pH and a resolution of 0.01 pH, very precise measurements can be taken with the school environmental meter. Because of the small and compact design, the school environmental meter is ideal for mobile use.

- Easy operation
- Accuracy of  $\pm 0.1$  pH
- Ideal for school facilities
- Removable electronics
- Small and handy meter
- Resolution of 0.01 pH

## **Specifications:**

## **Delivery scope:**

x Environmental Meter PCE-PH 15
x Screwdriver
x 1.5 V LR44 AG-13 batteries
x User manual