

KERN OBF133T241

The digital compound microscope for every laboratory with fixed, pre-centred Köhler illumination



Power supply

Input voltage: 100 - 240 V

Category

Brand: Optics

Packaging & shipping

Delivery: 24 h

Gross weight: 7,278 kg

Net weight: 9,050 kg

Shipping method: Parcel service

Construction

Dimensions (WxDxH): 394 x 185 x 377 mm

Eye distance (max): 75

Eye distance (min): 50

Optical system: Finite

Rotation of tube 360°: yes

Tube: Trinocular

Tube angle: 30°