



YUYANG INDUSTRIAL CO., LIMITED

China Manufacturer of Fire Testing Equipment

ASTM Flexure Test Fixture ISO8124-1 Aluminum Alloy Toys Safety Test Equipment





- **Product Details:**

- Place of Origin: **China**
- Brand Name: **YUYANG**
- Certification: **16 CFR 1500.48 16 CFR 1500.49**

- Model Number: YY1136
- **Payment & Shipping Terms:**
- Minimum Order Quantity: **1 set**
- Price: **Negotiation**
- Packaging Details: **Plywood Box**
- Delivery Time: **3-5 work days**
- Payment Terms: **T/T L/C Western Union**
- Supply Ability: **10 sets per month**
- Share to :

10 sets

ASTM Flexure Test Fixture ,ISO8124-1 Aluminum Alloy Toys Safety Test Equipment

Description:

To determine whether metallic wires or materials used for stiffening or retention of form in toys can be easily bent or broken. Potential hazards of sharp points & sharp edges.

Used to test the toys from the increase in rigidity or fixed shape of the metal wire or other metal materials used in the bending (or abuse) the process, whether due to fracture a risk of a sharp tip or sharp edges protruding. Abuse testing program is a toy. Made of stainless steel material. United States, Europe and China national standard test board size requirements vary.

Standards:

16 CFR 1500.48 16 CFR 1500.49

