

Description

According to Wolff Wilborn standard, this tester is designed to measure the coating surface hardness scrape-resistance, applied to experiment analysis or control production quality, suitable products including oil painting, furniture, car, cell phone, notebook, language translator, digital camera, etc. It is fast and easy to test single or double layer baking varnish.

Key words: abrasion testing machine, abrasion tester, abrasion test, abrasion, surface abrasion test, pencil hardness, JIS K5400 $_{\circ}$



Specifications

Pencil test angle	45°C angle for pencil and flat surface.
Pen nib load adjustment	0N ~ 10N
Test speed	10mm/sec (600mm/min)
Pensile specification	DIN standard test 6B ~ 9H (Optional)
	JIS standard test 6B ~ 9H (Optional)
Power	AC100 ~ AC240
Dimension	280 x 130 x 226 mm
Weight	4.8 kg

*Cometech reserves ultimate right to modify specification of the machine.



COMETECH TESTING MACHINES CO.,LTD. No.49,Ln,39,Lianucn Rd, Fengyuan Dist, Taichung City420,Taiwan Tel:886-4-25368881 / Fax:886-4-25368883 E-mail: sales@come-tech.com.tw URL: http://www.come-tech.com.tw

